



Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.



TRBC1/T-Cell Lymphoma Companion Panel

Methodology

Flow Cytometry

Test Description

Available as global and tech-only. Available as an add-on to the Standard Leukemia/Lymphoma Panel, and as stand-alone test (as described here) only for Breast capsules and fluids. Markers are CD3, CD4, CD7, CD8, CD25, CD26, CD30, CD45, CD279, and TRBC1 (10 markers).

Clinical Significance

This panel assesses T-cells for the presence of targetable antigens to guide therapeutic decisions. Expression of some markers may also help with differential diagnosis of various T-cell lymphomas or leukemias (including Mycosis fungoides/Sézary syndrome, anaplastic large cell lymphoma, angioimmunoblastic T cell lymphoma, etc.).

Specimen Requirements

- **Bone Marrow Aspirate:** 1-2 mL EDTA. Sodium heparin is acceptable. Lithium heparin or ACD (pale yellow/no gel separator) is not acceptable. Please provide recent CBC report.
- **Peripheral Blood:** 1-2 mL EDTA. Sodium heparin is acceptable. Lithium heparin or ACD (pale yellow/no gel separator) is not acceptable. Please provide recent CBC report.
- **Fresh Bone Marrow Core Biopsy:** 1-2cm core (length) tissue in RPMI
- **Fresh/Unfixed Tissue:** 0.2 cm³ minimum in RPMI
- **Fluids and FNAs:** Equal parts RPMI and specimen volume
- **CSF:** 1-2 mL recommended
- **NY Clients:** Please provide Date and Time of Collection.
- **Note:** Please exclude biopsy needles, blades, and other foreign objects from transport tubes. These can compromise specimen viability and yield, and create hazards for employees.

Storage & Transportation

Specimens should be received at NeoGenomics within 72 hours from collection to assure sample integrity and acceptable cell viability. **Note:** New York State samples must be received within 48 hours from collection per NYS requirements. Refrigerate specimen. Do not freeze. Use cold pack for transport, making sure cold pack is not in direct contact with specimen.

CPT Code(s)*

88184x1, 88185x9. Add 88188x1 for global

New York Approved

Yes

Level of Service

Technical, Global

Turnaround Time

24 hours

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way
Fort Myers, FL 33912
Phone: 239.768.0600/ Fax: 239.690.4237
neogenomics.com
© 2024 NeoGenomics Laboratories, Inc. All Rights Reserved.
All other trademarks are the property of their respective owners
Rev. 112424