

Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.





Non-Hodgkin's Lymphoma (NHL) FISH Panel

Alternative Name

NHL FISH

Methodology

FISH

Test Description

Probes: ALK (2p23) | BCL6 (3q27) | MYC (8q24) | CCND1/IgH t(11;14) | IgH (14q32) | IgH/BCL2 t(14;18) | MALT1 (18q21)

Probes may be ordered separately.

Disease(s): NHL

Clinical Significance

The NHL FISH panel is used for the detection of chromosome aberrations observed in lymphoma.

Specimen Requirements

- Bone Marrow Aspirate: 1-2mL Sodium Heparin Tube. EDTA tube is acceptable
- Peripheral Blood: 2-5mL Sodium Heparin Tube. EDTA tube is acceptable
- Fresh, Unfixed Tissue: Tissue in RPMI
- Bone Marrow/ Peripheral Blood Smear or Fresh Tissue Touch Preparation Slides: minimum 7 slides labeled with specimen type.
- Fluids: Equal parts RPMI to specimen volume.
- Fixed Cell Suspension: A client fixed cell suspension may be submitted for testing as long as it is received in 3:1
 Methanol:Glacial Acetic Acid.
- Paraffin Block: H&E slide (required) plus paraffin block. Circle H&E for tech-only.
- Cut Slide: H&E slide (required) plus 12 unstained slides cut at 4 microns for panel, or 2 unstained slides cut at 4 microns for single abnormalities. Circle H&E for tech-only.
- **Note:** Please exclude biopsy needles, blades, and other foreign objects from transport tubes. These can compromise specimen viability and yield, and create hazards for employees.

Storage & Transportation

Refrigerate specimen. Do not freeze. Use cold pack for transport, making sure cold pack is not in direct contact with specimen. For fresh samples: ship same day as drawn whenever possible; specimens <72 hours old preferred.

CPT Code(s)*

88374x7 automated. Codes may differ if manual analysis is performed.

New York Approved

Yes

Level of Service

Technical, Global

Turnaround Time 3-5 days	
*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole	

responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way Fort Myers, FL 33912

Phone: 239.768.0600/ Fax: 239.690.4237

neogenomics.com

 $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \ \,$ $\ \$ $\$ $\$ $\ \$ $\$ $\$ $\ \$ $\$ $\$ $\ \$ $\$ $\$ $\$ $\$ $\ \$ $\$ $\$ $\$ $\$ $\ \$ $\$ $\$ $\$ $\$ $\ \$ $\$

Rev. 041824