

# Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.





# Extended Leukemia/Lymphoma Panel - 31 markers

#### Methodology

Flow Cytometry

## **Test Description**

Available as global and tech-only. Markers are CD2, CD3, CD4, CD5, CD7, CD8, CD10, CD11b, CD11c, CD13, CD14, CD15, CD16, CD19, CD20, CD23, CD33, CD34, CD38, CD41, CD45, CD56, CD64, CD71, CD117, CD138, CD235a, FMC-7, HLA-DR, kappa, and lambda.

### **Clinical Significance**

For diagnosis of leukemia, lymphoma, plasma cell neoplasms, and evaluation of myeloid maturation.

#### **Specimen Requirements**

- Bone Marrow Aspirate: 1-2 mL EDTA. Sodium heparin is acceptable. Lithium heparin or ACD (pale yellow/no gel separator) is not acceptable. Please provide recent CBC report.
- **Peripheral Blood:** 1-2 mL EDTA. Sodium heparin is acceptable. Lithium heparin or ACD (pale yellow/no gel separator) is not acceptable. Please provide recent CBC report.
- Fresh Bone Marrow Core Biopsy: 1-2cm core (length) tissue in RPMI
- Fresh/Unfixed Tissue: 0.2 cm3 minimum in RPMI
- Fluids and FNAs: Equal parts RPMI and specimen volume
- CSF: 1-2 mL recommended
- NY Clients: Please provide Date and Time of Collection.
- **Note:** Please exclude biopsy needles, blades, and other foreign objects from transport tubes. These can compromise specimen viability and yield, and create hazards for employees.

#### Storage & Transportation

Specimens should be received at NeoGenomics within 72 hours from collection to assure sample integrity and acceptable cell viability. Note: New York State samples must be received within 48 hours from collection per NYS requirements. Refrigerate specimen. Do not freeze. Use cold pack for transport, making sure cold pack is not in direct contact with specimen.

#### CPT Code(s)\*

88184(x1), 88185(x30). Add 88189(x1) for global.


#### **New York Approved**

Yes

#### **Level of Service**

Technical, Global

#### **Turnaround Time**



Please direct any questions regarding coding to the payor being billed.

<sup>\*</sup>The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

\*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way Fort Myers, FL 33912

Phone: 239.768.0600/ Fax: 239.690.4237

neogenomics.com

Rev. 042524