

Background

Gene fusions are major drivers of cancer and their accurate detection is key for supporting diagnosis and therapy selection. There are more than 65,000 annotated gene fusion events. Current fusions detection strategies are limited to amplicon panels querying only a small number of genes, or whole exome capture based assays, which reduces sensitivity. Hence, multiple small panels are required to capture all relevant fusions, which is prohibitive for small biopsies like fine needle aspirates. Hence, a clinical grade assay able to accurately detect a comprehensive set of most clinically relevant gene fusions for solid tumors is needed.

Methods

We isolated total nucleic acid (TNA) from 141 FFPE tumor specimens or a control sample (Sequencing) and performed pair-ended, strand-specific hybridization-based RNA sequencing on Next-Seq 500 or Novaseq 6000 platforms. We selected 252 fusion genes from NCCN and WHO guidelines, published clinical studies, and the 120 most frequent curated fusions in solid tumors from COSMIC (v91). These fusions are clinically relevant to most frequent cancers, including breast, colorectal, lung, lymphoma, pancreatic, prostate, salivary gland, sarcomas, and thyroid cancers. Chimeric probes were synthesized targeting fusion RNA sequences for 2230 selected breakpoints and exon junction regions of aberrant RNAs, including EGFRvIII, MET exon 14 skipping, ARV7 and ARV9. We also targeted the full coding region of 27 genes whose change in expression may suggest translocations and mutations, or have diagnostic value. Fusions were called by a custom pipeline using 3 fusion callers and a machine-learning algorithm. The PCR-based Archer FusionPlex assays and RT-PCR followed by Sanger were used as orthogonal validation methods.

Results

- This single RNA fusion assay can detect at least 1194 unique fusions pairs involving 1104 genes with 250 core fusion genes compared to 63 or 47 fusion genes from other commercially available NGS assays.
- CLIA assay validation was conducted on 146 unique clinical FFPE samples from various tumor types.
- 206 unique genes were detected in 1 or more gene fusions. The assay detected 116/121 high confidence fusions reported in our CLIA Archer Sarcoma and Comprehensive Thyroid & Lung (CTL) assays. 41 additional fusions not detected by Archer were confirmed by Sanger sequencing.
- The assay has 95.8% sensitivity (141/147) and 100% specificity, as all fusions were confirmed by either orthogonal assay.
- We also detected MET exon 14 skipping, EGFRvIII with 100% sensitivity and 100% specificity. Importantly, in two samples we detected MET exon 14 skipping not predicted from DNA mutation analysis showing the sensitivity of the approach.
- 55 novel fusions were detected by capturing only one of the partner genes.

Conclusions and Future Directions

We developed a novel and efficient breakpoint targeted fusion detection RNA-seq assay from extracted TNA from FFPE samples that can comprehensively profile thousands of clinically actionable RNA fusions and aberrant RNAs in solid tumors. Future directions include validating the built-in gene expression module to complement even further the characterization of solid tumors by integrating the full spectrum of RNA alterations from this assay with DNA sequencing assays reporting mutations/CNVs thus enhancing Precision Oncology alternatives for cancer patients.

Data

Gene Fusions Detection: RNA vs. DNA Detection




Figure 1: Schematic representation of gene fusions breakpoints and the advantage for fusion detection on mRNA.

Solid Tumor Relevant Gene Fusions and Aberrant Transcripts

250 clinically relevant fusion genes											
ABL1	C11orf95	CREB1	ESRP1	HERPUD1	MAST2	NOTCH2	PML	SET	TFE3		
ACSL3	CAMTA1	CREB3L1	ETV1	HEY1	MEAF6	NPM1	POU5F1	SLC34A2	TFG		
ACTB	CANT1	CREB3L2	ETV4	HIP1	MET	NRA43	PPARG	SLC45A3	THADA		
ACTL6A	CAPZA2	CREBBP	ETV5	HMGGA2	MET e14*	NRG1	PRCC	SN01	THRAP3		
AFDN	CARS1	CRTC1	ETV6	HMGN2P46	MKRN1	NTRK1	PRKACA	SNURF	TMPRSS2		
AFF1	CBFA2T3	CRTC3	EWSR1	HNRNP2A81	MLL1	NTRK2	PRKAR1A	SRF	TP63		
AFF3	CCDC170	CSF1	EZR	IL2RB	MLL10	NTRK3	PRKD1	SRGAP3	TPM3		
AFF4	CCDC6	CTLA4	FAM131B	IRF4	MLL11	NUP214	PRKD2	SS18	TPM4		
AHRR	CENB3	CTNNA1	FEV	ITK	MLL3	NUP98	PRK3	SSX1	TPR		
AKAP9	CCND1	CTNNA1	FGFR1	JAK2	MN1	NUTM1	PTPRK	SSX2	TRIM24		
AKT3	CCND2	DDIT3	FGFR2	JAZF1	MPRIP	NUTM2A	RAD51B	SSX4B	UBTF		
ALK	CCND3	DDX3X	FGFR3	KAT6A	MRTFB	NUTM2B	RAF1	STAT6	USP6		
AR	CD274	DDX5	FGFR4	KDM5A	MSH2	OMD	RANBP2	STIL	VTG1A		
ARV7*	CD74	DHH	FLI1	KIAA1549	MYB	PAN3	RARA	STRN	WDFY2		
ARV9*	CDH11	DNAJB1	FLT3	KIF5B	MYBL1	PATZ1	RASGEF1A	SUZ12	WIF1		
ARID1A	CDK4	DUX4	FOXO1	KIT	MYC	PAX3	RET	SYK	WT1		
ASPSR1	CDK6	EGFR	FOXP1	KLK2	MYH9	PAX7	RHEBL1	TACC3	WWTR1		
ATF1	CDKN2D	EGFRvIII*	FRK	KMT2A	MYLK	PAX8	ROS1	TAF15	YAP1		
ATIC	CHCHD7	ELK4	FUS	KNL1	NAB2	PBX1	RP56KC1	TAL1	YWHAE		
AXL	CIC	ELL	GLI1	KRAS	NCOA1	PCM1	RSP03	TBL1XR1	ZNF444		
BCOR	CIITA	EML4	GLIS1	LIFR	NCOA2	PDGFBR	RUNX1	TCF12			
BCR	CLTC	EPC1	GLIS2	LMNA	NCOA4	PDGFRA	RUNX1T1	TCF3			
BRAF	CNBP	EPS15	GLIS3	LPP	NDRG1	PDGFRB	SDCA	TCF7L2			
BRCA1	COA5	ERBB2	GNAS	MAML1	NFATC2	PHE1	SEC31A	TEAD1			
BRCA2	COL1A1	ERG	GOPC	MAML2	NFB	PIK3CA	SEPTIN6	TEAD2			
BRD4	COL1A2	ESR1	HAS2	MAST1	NOTCH1	PLAG1	SEPTIN9	TEAD3			

Table 1. the actionable genes reported on the universal fusion comprehensive fusions panel.

RNA Fusions Capture-based Targeted RNA-seq Assay Workflow




Figure 2: A new Solid Tumor Gene Fusion RNA-Seq assay workflow. Starting with FFPE total nucleic acid extraction, the assay involves an RNA fusion sequence- targeted, hybridization capture-based, stranded pair-ended RNA-sequencing. Key step resides on custom based design of capture baits mapping to the chimeric RNA sequence as deduced from chromosomal breakpoints in annotated private and public databases.

CLIA Validation Samples Cohort




Figure 3: Disease representation of samples analyzed in the CLIA/CAP assay validation.

CLIA Validation Fusion Genes




Figure 4: Fusion genes frequency. From 206 different genes, partners on the fusions detected in this validation, 51 genes were found as partners of >1 unique gene fusions.

New Fusions Confirmation: RT-PCR +Sanger Sequencing




Figure 5: Novel fusions confirmation. All novel fusions detected based on the capture by a probe hybridizing on a known fusion gene were validated as the fusion shown here.

Assay Performance

Variant Type	Tested tumor cellularity	Sensitivity	Specificity
All RNA Variants	20%	95.8%	100.0%
RNA Fusions (gene A -gene B)	20%	95.5%	100.0%
Splicing/mutant transcripts (MET exon 14 skipping)	20%	100.0%	100.0%
Mutant derived transcripts (EGFRvIII)	20%	100.0%	100.0%

Table 3. Performance results from accuracy experiments.

RNA-seq vs. DNA-seq Mutations Testing for Splicing Alterations





Figure 6: Diagram of splicing inducing MET DNA mutations and clinical results from RNA-seq vs. DNA sequencing results for 2 MET exon 14 Positive samples also analyzed by a CLIA validated NGS DNA test. 83 total samples were tested for MET exon 14 skipping.

New Fusions and New Partner Fusion Genes Detection



- Known Fusions: 2230 custom capture probes for known gene fusions
- Novel fusions: one sided capture of novel fusion partners
- Probes for breakpoints on 1104 unique partner genes enabling the capture of fusions to anywhere on the transcriptome.

Figure 7: Coverage by capture probes on new fusions confirmed by RT-PCR +Sanger sequencing. Percentage of fusions targeted by the probes at both sides or one side of the fusion is indicated.