

## Introduction

In recent years, novel immunotherapy strategies have shown remarkable results in treating various advanced cancers. However, the efficacy of immunotherapeutics has been dictated by the tumor microenvironment (TME), immune tumor infiltrates, immune checkpoint protein expression, and molecular tumor profiles leading to widely variable clinical outcomes and responses in patients. The dynamic and spatially heterogeneous TME is made up of malignant and nonmalignant cells such as endothelial cells, cancer-associated fibroblasts, immune cells, adipose cells, and neuroendocrine cells in addition to vascular and lymphatic networks and the extracellular matrix. The local milieu as well as complex tumor-TME interactions vary widely between tumor of the same type and also within the same tumor and therefore highlights the need for enhanced analysis of tissue-immune content and their spatial context for predicting response and biomarker discovery.


In this study, we applied a multi-faceted, highly multiplexed tissue analysis approach to quantitate and better characterize the spatial arrangement of key immuno oncology protein markers in a pan-cancer cohort of up to 35 FFPE samples from breast, head and neck, prostate, non-small cell lung cancer (NSCLC), endometrial and colorectal indications using NanoString GeoMx™ Digital Spatial Profiler (DSP) and MultiOmyx™ Hyperplexed Immunofluorescence (IF) assay.

For DSP analysis, each sample was profiled for up to 52 protein biomarkers included on NanoString Human Immuno oncology (IO) protein panels. Selection of regions of interest will be guided by H&E staining, fluorescent markers (CD45, PanCk and Syto 13). Protein profiling of tumor and TME regions was achieved by segmenting by PanCk+/PanCk- following by collection of indexed oligonucleotides and digital counting using NanoString nCounter system. Data normalization using endogenous controls across the ROIs were performed to understand the spatial distribution of these IO biomarkers within the TME.

An adjacent section of each tumor sample was stained and analyzed using MultiOmyx technology. MultiOmyx is a proprietary IF platform for the visualization and characterization of up to 60 protein biomarkers in a single FFPE section. Using the MultiOmyx IF assay in combination with proprietary algorithms, we studied the expression and spatial distribution of 14 immuno oncology (IO) markers to characterize different subtypes of immune cells and immune checkpoint inhibitors within the TME in the FFPE tumor samples.

Combination of MultiOmyx IF with NanoString DSP provides a complementary and powerful solution to study the IO markers within the TME. Using the DSP nCounter counts along with MultiOmyx cell classification and spatial analysis results, the researchers can simultaneously profile the high-plex protein expression, characterize the immunophenotypes and visualize the spatial distribution of tumor infiltrating immune cells at single-cell resolution within the TME.

## Overview of Study Workflow




## Key Findings

In this study, NanoString technologies GeoMx DSP and MultiOmyx, both platforms offered by NeoGenomics Laboratories were utilized to perform expression profiling and immune cells phenotyping in breast, head and neck, prostate, lung, endometrial and colorectal cancer samples.


- Correlation coefficient was calculated between GeoMx DSP counts and the positive cell densities measured by MultiOmyx IF assay for the 10 markers that showed enough positivity in this study (>100 positive cells per ROI). Direct correlation was observed for 8 out of 10 markers.
- This integrated analysis by both technologies provides comprehensive understanding of the immune landscape in oncology FFPE tissues.
  - GeoMx DSP provides a high multiplexing, non-destructive, spatial approach to overall expression profiling of the 52 IO markers from FFPE samples tested in the study.
  - MultiOmyx is a robust approach for immunophenotyping and characterizing the co-expression/co-localizing of the biomarkers at single cell level.

## Expression Profiling using NanoString GeoMx DSP




**Figure 2. Analysis of biomarker profiling using DSP in FFPE samples.** A. A representative DSP image showing ROI selection and segmentation using PanCk. B. Normalized counts of PanCk and CD3 in PanCk+ (Tumor) and PanCk- (TME) regions. C. Reproducibility of DSP biomarker expression between when tested using different analysts. D and E. Repeatability of DSP biomarker expression between ROIs within a sample (D) and between serial sections of the sample (E). F. Heatmap of the biomarker density results (for a subset of biomarkers) in the all tissue areas.

## Immunoprofiling using MultiOmyx Immunofluorescence Assay



**Figure 3. Immunoprofiling of the FFPE Tissues by MultiOmyx.** A. Quantification and classification results of the immune phenotypes in the study. B. Heatmap of the biomarker density results in the all tissue areas. C. Heatmap of the co-expressions in the all tissue areas. D. Representative color overlay images from the sample S2000114-SEQ-00007. This sample showed high expression of immune markers by both MultiOmyx and DSP.

## Correlation Study



**Figure 4. Correlation Between NanoString DSP and MultiOmyx IF.** A. Scatter plots showed the correlation between NanoString DSP counts and the densities of positive cells measured by MultiOmyx IF assay. B. Correlation coefficient was calculated for the 10 markers with enough positivity (>100 positive cells per ROI) in this study. Direct correlation was observed for 8 markers.