

# Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.





# **NeoTYPE® Liver/Biliary Tumor Profile**

#### **Alternative Name**

Liver/Biliary Tumor Profile

### Methodology

Molecular

# **Test Description**

The NeoTYPE Liver/Biliary Tumor Profile analyzes 38 biomarkers through a combination of next-generation sequencing (NGS), FISH, and IHC as listed below. Test orders include summary interpretation of all results to help guide treatment decisions. If Pan-TRK IHC is expressed or equivocal, reflex to either NTRK NGS Fusion Panel (Default) or NTRK 1-3 FISH Panel will be added. A microsatellite instability (MSI) NGS result of "indeterminate" will create a reflex to MSI by PCR as long as the tumor percentage is ?40% and paired normal tissue is available.

- NGS (32 genes + 2 biomarkers): AKT1, APC, ARID1A, ATM, BAP1, BRAF, CDKN2A, CTNNB1, EGFR, ERBB2, ERBB4, FGFR1 (mutations), FGFR2 (mutations), FGFR3 (mutations), HRAS, IDH1, IDH2, KRAS, MET, Microsatellite Instability (MSI), NF1, NOTCH1, NRAS, PBRM1, PIK3CA, PTEN, SMAD4, SMO, SRC, TERT Promoter, TP53, TSC1, TSC2, Tumor Mutation Burden (TMB)
- FISH (3 FISH): FGFR2 rearrangement, MET, PTEN (tech-only available)
- IHC (2 biomarkers): PD-L1 LDT, Pan-TRK (tech-only available for PD-L1)

# Clinical Significance

The NeoTYPE Liver/Biliary Tumor Profile characterizes liver/biliary tumors of any histological subtype for the most significant genetic changes relevant to therapy decisions, prognosis, and clinical research. It is appropriate for patients with newly-diagnosed or recurrent disease, and for patients with resistant disease to explore options in clinical trials.

#### **Specimen Requirements**

• FFPE tissue: Paraffin block preferred. Please use 10% buffered formalin fixative. Do not use zinc fixatives.

#### **Storage & Transportation**

Use cold pack for transport, making sure cold pack is not in direct contact with specimen.

## CPT Code(s)\*

81445x1, 88377x3, 88360x1, 88342x1; add 81479x1 if reflexed to NTRK NGS Fusion Panel (default) or 88374x3 automated (88377x3 manual) if reflexed to NTRK 1-3 FISH Panel

#### Medicare MoIDX CPT Code(s)\*

81479

# **New York Approved**

Yes

# **Level of Service**

Global

# **Turnaround Time**

14 days

Please direct any questions regarding coding to the payor being billed.

<sup>\*</sup>The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

\*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way Fort Myers, FL 33912

Phone: 239.768.0600/ Fax: 239.690.4237

neogenomics.com

© 2024 NeoGenomics Laboratories, Inc. All Rights Reserved. All other trademarks are the property of their respective owners

Rev. 052024