

Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.





EBER

Alternative Name

Epstein-Barr Virus by ISH

Methodology

In Situ Hybridization (ISH)

Test Description

This probe set labels all latent EBV-infected cells, including EBV-positive lymphoblastoid cell lines and EBV infected B-cell immunoblasts in infectious mononucleosis. It also reacts with EBV-associated undifferentiated nasopharyngeal carcinomas and with Reed-Sternberg cells in almost all EBV-associated Hodgkin lymphoma cases. Global interpretation is available on head and neck specimens only; tech-only testing is available for all samples.

Clinical Significance

The Epstein-Barr virus (EBV) probe demonstrates latent EBV infection by hybridizing to abundantly expressed EBER transcripts which are concentrated in the nuclei of latently infected cells.

Specimen Requirements

- A formalin-fixed, paraffin-embedded (FFPE) tissue block is preferred specimen type or
- One (1) unbaked, unstained slide for H&E staining (required) and three (3) positively charged unstained slides, all cut
 at 4-5 microns
- Block and slide identifiers should be clearly written and match exactly with the specimen ID and specimen labeling as noted on the requisition.

Storage & Transportation

Use cold pack for transport. Make sure cold pack is not in direct contact with specimen.

CPT Code(s)*

88365 x 1 or 88364 x 1

New York Approved

Yes

Level of Service

Stain Only, Global

Turnaround Time

Global (head and neck only): 3 days, Tech-Only (stain only): 48 hours

Please direct any questions regarding coding to the payor being billed.

^{*}The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way Fort Myers, FL 33912

Phone: 239.768.0600/ Fax: 239.690.4237

neogenomics.com

Rev. 042024